

# K5 cryolaser®



## K5 cryolaser®

Brown skin spots can also be lightened by laser application, but removal of brown spots by the laser is sometimes incomplete and a brown circular spot all around the area of the lightened spot can remain.



**K5 cryolaser®** is an anti-brown spots skin lipogel **expecially formulated** for a use combined with the laser application.

The aim of its use is the **lightening of the thin brown spot remaining after laser application.**

**K5 cryolaser® does not contain salicylic acid.**

After laser application the skin is damaged by it and the use of salicylic acid is not advisable and not necessary.